Vitalograph micro Regulatory Notices

Warning and Advisory Notices:

Note: Please read all this information before using this Vitalograph device.

A full set of instructions, including cleaning instructions, is available at www.vitalograph.co.uk.

- The Vitalograph micro in association with Vitalograph Reports software is a spirometer designed for lung function testing in a variety of environments, e.g. hospital wards, occupational health centres and private homes.
 - Safety or Bacterial Viral Filter disposable mouthpieces must be used. The use of a disposable noseclip is also recommended.
 - Take care not to block the mouthpiece with the tongue or teeth. A 'spitting' action or coughing will give false readings.
 - Time zero is determined using the back-extrapolated method, from the steepest part of the curve.
 - An accuracy check should be performed daily with a precision syringe as recommended by the ATS and ERS. It is recommended to perform an accuracy check after annual maintenance, after cleaning or disinfecting the device, after adjusting calibration or if the device or flowhead has been dropped.
 - If the device is used longer than its specified life, the accuracy of the device may deteriorate.
 - Symptoms must take precedence over device measurements.
 - Store in a clean dry place.
 - Cleaning and disinfection:
 - The outer surfaces should be cleaned once a week, more often if necessary. The use of an ordinary alcohol wipe is recommended, with special attention to the mouthpiece area.
 - The flowhead may be disinfected using sodium dichloroisocyanurate solution.
- Refer to User Training Manual for full cleaning instructions.
-  The device must be taken to separate collection at the product end-of-life.
-  Do not dispose of these products as unsorted municipal waste.

Warranty

Your Vitalograph device is guaranteed for one year*. Replace it if it is faulty, otherwise replace the unit every five years.

Any unauthorised changes to Vitalograph device hardware or software may compromise product safety and/or data and as such Vitalograph cannot be held responsible and the device will no longer be supported

*Excepting accidental / transit damage or inappropriate use of the device.

FDA Notice:


Caution: Federal Law restricts this device to sale by, or on the order of a physician.

CE Notice

Marking by the symbol  indicates compliance of the Vitalograph micro to the Medical Devices Directive of the European Community. Such marking is indicative that the Vitalograph micro meets or exceeds the referenced technical standards.

Medical Devices may be affected by cellular telephones and other personal or household devices not intended for medical facilities. It is recommended that all equipment used near the Vitalograph product comply with the medical electromagnetic compatibility standard and to check before use that no interference is evident or possible. If interference is suspected or possible, switching off the offending device is the normal solution, as is required in aircraft and medical facilities.

The Vitalograph device is intended for use in the electromagnetic environment specified below. The customer or user of the system should assure that it is used in such a system.

- Floors should be wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be less than 30%.
- Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
- Interference may occur in the vicinity of equipment marked with the following symbol- 

Recommended separation distances between portable and mobile RF communication equipment and the Vitalograph micro

The Model 6300 micro is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the micro can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the micro as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz d = 1.2√P	80 MHz to 800 MHz d = 1.2√P	800 MHz to 2.5 GHz d = 2.3√P
0.01	0.1m	0.1m	0.2m
0.1	0.4m	0.4m	0.7m
1	1.2m	1.2m	2.3m
10	3.7m	3.7m	7.4m
100	11.7m	11.7m	23.3m

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Technical Specification

Flow detection principle	Fleisch type pneumotachograph
Back pressure	Less than 0.1kPa/L/second @ 14L/s, complies with ISO 26782
Volume detection	Flow integration sampling @ 100Hz
Maximum test duration	90 seconds
Maximum displayed volume	10L
Volume accuracy	better than ±3% (or 0.05L) Max. Volume 9.99L Min. Volume 0.01L
Flow accuracy when operated in operating temperature range conditions	Flow ±10% Max. flow rate ±16 L/s Min. flow rate ±0.02 L/s
Linearity	Better than ±3%
Voltage	5V USB /6V Battery (4 x 1.5V Alkaline AAA)
Operating temperature range	17–37°C (ISO26782 limits), Design limits 10–40°C
Performance standards the Vitalograph micro meets or exceeds	EN ISO 23747:2009, ISO 26782:2009, ATS/ERS:2005
Safety standards	EN 60601 (IEC 60601)
Medical safety standards	Medical Devices Directive 2007/47/EC
QA/GMP standards	ISO 13485:2003, FDA 21 CFR 820, CMDR SOR 98/282 & JPAL
Storage temperature	0–50°C
Storage relative humidity	10%–95%
Communications	USB

Declaration of Conformity

Product: Model 6300 Vitalograph micro

Vitalograph hereby ensures and declares that the above product associated with this user manual, is designed and manufactured in accordance with the following QMS regulations and standards:

- European Medical Devices Directive (MDD) 93/42/EEC, as amended.
- This device is classified as 2a per Annex IX of the MDD also meets the provisions of the Essential Requirements, Annex I, via compliance with Annex II of the Medical Devices Directive as per Article 11, section 3a, excluding point 4 of Annex II.
- Canadian Medical Device Regulation (CMDR SOR/98-282)
- FDA Quality System Regulation (QSR) 21 CFR 820.
- EN ISO 13485:2003. Medical devices. Quality management systems. Requirements for regulatory purposes.

Certifying Body: British Standards Institute (BSI).
(For 93/42/EEC and CMDR).

BSI Notified Body #: 0086

Certificate Nos. CE 00772, CE 85553, MD 82182, FM 83550

Signed on behalf of Vitalograph (Ireland) Ltd.

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Firmado en nombre de Vitalograph (Ireland) Ltd.