

This is a Customer Header for Company XYZ, Dept ABC

### Subject Information

Subject Name: John Doe	Subject ID: 564565
Gender: Male	Age: 37
Height: 175 cm	Regression Set: Temp
Weight: 70 kg	BMI: 22.9


### Test Information

Device ID: 1234567890	Report Date: 09-Jul-2008 13:28
Device Firmware Ref: 04108	Device Software Ref: V1.06
Number of Blows: 10	Number of Good Blows: 5
Best FEV1 within: 0.04 L	Best FEV6 within: 0.2 L

### Results


Parameter	Normal Predicted	Test 1	Test 2	Test 3	Best Test	% Pred
FEV1 (L)	<b>3.96</b>	2.25	2.21	2.16	<b>2.25</b>	<b>57</b>
FEV6 (L)	<b>4.78</b>	2.48	2.28	2.16	<b>2.48</b>	<b>52</b>
FEV1/FEV6 (ratio)	<b>0.82</b>	0.90	0.96	0.99	<b>1.00</b>	<b>110</b>

### Obstructive Index



FEV1 Measured/FEV1 Predicted

### COPD (GOLD) Classification




### Interpretation



Lung Age: 89 years

Obstructive Index: 57%

Interpretation: Normal, not COPD  
Possible Restriction



### Comments

Test Report

**Signature:** \_\_\_\_\_

**Date:** \_\_\_\_\_