

Clinic & Device Information

Device ID:	0000110264	Device Software Ref:	40915 V1.00
Clinic ID:	123	Clinic:	Paris Clinic


Subject Information

Subject Name:	Peter B Wang	Subject ID:	N1000456
Gender:	Male	Age:	46
Height:	170 cm	Date of Birth:	01/03/1963
Weight:	100 kg	BMI:	34.6
Set FEV1:	2.50(L)		

Results Summary

Highest FEV1:	2.64 (L)	Start Date:	20/07/2009 08:57:52
Lowest FEV1:	1.66 (L)	End Date:	23/07/2009 09:39:59
Mean FEV1:	2.26 (L)	Number of Days:	4
Highest FEV6:	3.11 (L)	Number of Sessions:	15
Lowest FEV6:	1.95 (L)	Good Sessions:	14
Mean FEV6:	2.64 (L)	Bad Sessions:	1
Highest FEF 25-75:	3.67 (L/sec)		
Lowest FEF 25-75:	1.52 (L/sec)		
Mean FEF 25-75:	2.40 (L/sec)		

FEV1 Chart



Comments

Chronic lung disease

Signature: _____

Date: _____



Lung Monitor

Your header Respiratory Partner

FEV1 Zones

[0] Green	[>2.38]
[1] Yellow	[2.13-2.38]
[3] Red	[<2.13]

Detailed Data Table

#	Date/Time	Session Best FEV1 (L)	FEV1 Zone	Session Best FEV6 (L)	Session Best FEF 25-75 (L/sec)	Blows	Session Quality
1	20/07/2009 08:57:52	2.53	0	2.83	3.67	3	Good
2	20/07/2009 09:07:17	2.43	0	2.90	2.43	2	Good
3	20/07/2009 09:12:43	2.62	0	2.99	2.77	3	Good
4	20/07/2009 09:15:13	2.64	0	3.06	2.72	3	Good
5	20/07/2009 09:48:34	2.64	0	3.11	2.60	2	Good
6	20/07/2009 09:50:34	2.64	0	3.07	2.69	3	Good
7	20/07/2009 12:56:51	2.48	0	2.76	2.79	3	Good
8	20/07/2009 19:00:27	1.99	3	2.01	2.40	1	Bad
9	20/07/2009 20:54:01	2.31	1	2.54	2.46	1	Good
10	21/07/2009 07:31:26	1.69	3	2.02	1.68	3	Good
11	21/07/2009 20:46:27	2.07	3	2.68	1.81	3	Good
12	22/07/2009 07:18:48	1.66	3	2.25	1.52	3	Good
13	22/07/2009 22:03:44	2.23	1	2.95	2.09	3	Good
14	23/07/2009 07:33:55	1.78	3	1.95	1.85	2	Good
15	23/07/2009 09:39:59	2.23	1	2.46	2.52	3	Good